



Hingeless ankle

Code	Material	Screw	Weight (max), kg	Max patient's weight, kg
S8-M10	Stainless steel	Stainless steel	0.13	125
S8T-M10	Titan	Stainless steel	0.10	125

Ankle with ankle joint hinge

Code	Material	Dimensions: length/width/height (max), mm	Weight (max), kg	Max patient's weight, kg
S1025M1	Stainless steel	90/45/75	0.50	125
S1029M1	Stainless steel	100/45/80	0.45	125

Lower leg holding module

Code	Material	Dimensions: length/width/height (max), mm	Weight (max), kg	Max patient's weight, kg
S3-200	Stainless steel/aluminum	40/200	0.20	100
S3T-200	Titanium/aluminum	40/200	0.15	100

Upper leg holding module

Code	Material	Dimensions: length/width/height (max), mm	Weight (max), kg	Max patient's weight, kg
S3-400	Stainless steel/aluminum	40/400	0.3	100
S3T-400	Titanium/aluminum	40/400	0.26	100

Clamp adapter

Code	Material	Dimensions: length/width/height (max), mm	Weight (max), kg	Max patient's weight, kg
S 21	Stainless steel	40/400	0.11	125
S 52	Titanium	40/400	0.07	125

Double adapter

Code	Material	Dimensions: length/width/height (max), mm	Weight (max), kg	Max patient's weight, kg
S72-L32	Stainless steel	40/32	0.13	125
S72T-L32	Titanium	40/32	0.09	125
S72-L45	Stainless steel	40/45	0.15	125
S72T-L45	Titanium	40/45	0.1	125
S72-L60	Stainless steel	40/60	0.17	125
S72T-L60	Titanium	40/60	0.12	125
S72-L75	Stainless steel	40/75	0.19	125
S72T-L75	Titanium	40/75	0.13	125

Cased adapter

Code	Material	Dimensions: length/width/height (max), mm	Weight (max), kg	Max patient's weight, kg
S63	Stainless steel	70/33	0.11	125
S100	Titanium	63/31	0.06	125

Cased adapter

Code	Material	Dimensions: length/width/height (max), mm	Weight (max), kg	Max patient's weight, kg
S41	Stainless steel	117/85/61	0.17	125
S41T	Titanium	117/85/61	0.14	125

Cased adapter

Code	Material	Dimensions: length/width/height (max), mm	Weight (max), kg	Max patient's weight, kg
S89	Stainless steel	117/85/67	0.16	125
S89T	Titanium	117/85/67	0.14	125